

AET/APT

Tourniquets & Accessories



Designed to perform



Quality Assured

AET

Electronic Tourniquet

State of the art touch screen technology makes the AET so easy to use.



User friendly, easy to clean with low maintenance costs, the AET will be found ideal for Theatre, A & E and Pain Relief Clinics.

- Independently operated single or dual channel option colour coded



SAFETY FEATURES

Automatic timers, a screen safety lock and audible alarms provide complete safety when in use.

- State of the art touch screen technology
- Integral power supply
- Leak compensation facility
- Automatic self test routines
- User friendly and simple to use controls



Independently operated single and dual channel facility and an integral power supply.



- Automatic timers
- Screen safety lock
- Audible alarms
- Patient's data memory store
- Low maintenance
- Supplied with mobile stand and storage facilities



PRESSURE

Both channels can be used simultaneously, whilst remaining independently operated. Each channel has two secondary pressure gauges, one for pre-selection of cuff pressure, and the other showing actual cuff pressure on inflation.



APT

Tourniquet

Highly acclaimed for its reliability and user friendly design

The APT Tourniquet is a dual channel self-compensating tourniquet inflation unit. The separate channels are indicated as A & B and are colour coded blue and black. The colour coding is applied also to the cuff supply hoses to avoid confusion during use.



- Precision controlled quality pressure regulators
- Automatic timers independently activated through inflate/deflate mode
- Twin gauges ensure accurate pressure settings



Available with both piped service and 'E' size air cylinder options.

ACCESSORIES

RHYS DAVIES EXSANGUINATOR

The Rhys Davies Exsanguinator is a simple, but effective piece of equipment which allows the exsanguination of limbs for surgery without subjecting the limb to high pressures experienced with some other methods of exsanguination.



SERVICING

Understanding the demands required of the modern hospital environment, it is imperative these precision devices are maintained to the highest standards.



Supplied with a unique 3 year warranty option, Anetic Aid Skilled Personnel are constantly available for technical support and can offer routine maintenance at competitive rates.

ACCESSORIES

A complete selection of accessories are available including various size tourniquet cuffs, Y connectors and safety connectors.



New to our range is the Ultra Clean Tourniquet Cuffs, manufactured with welded seams, these totally sealed units are ideal for the modern hospital environment.



TOURNIQUET CUFF COVERS

Responding to MDA safety notice MDA SN 1999[26]: we are pleased to offer within our range the new "TOURNIQUET CUFF COVERS". From prepping, to completion of the procedure, this product helps prevent solution from reaching the tourniquet from either side, which eliminates skin irritation and allows a clean tourniquet for the next procedure.

ANETIC AID

The Theatre Equipment Specialists

AET/APT

Tourniquets & Accessories



ITEM

Cat No

AET ELECTRONIC TOURNIQUET

AET Electronic Tourniquet (c/w mobile stand & storage)

40050

APT TOURNIQUET CONTROL UNITS

APT MK3 Tourniquet Control Unit (piped air, 4 bar)

40000

APT MK3 Tourniquet Control Unit (piped air, 7 bar)

40010

APT MK3 Tourniquet Control Unit (cylinder model)

40020

MOBILE STANDS

Mobile Stand

40100

Mobile Stand (c/w yoke attachment, cylinder model)

40110

TOURNIQUET CUFFS

BLADDER SIZE

305 mm x 60 mm

Single

40200

305 mm x 120 mm

Double

40210

460 mm x 100 mm

Single

40220

460 mm x 140 mm

Double

40230

610 mm x 100 mm

Single

40240

610 mm x 150 mm

Double

40250

860 mm x 100 mm

Single

40260

1070 mm x 100 mm

Single

40270

TOURNIQUET CUFFS - ULTRA CLEAN

BLADDER SIZE

460 mm x 100 mm

Single

40220-U

610 mm x 100 mm

Single

40240-U

860 mm x 100 mm

Single

40260-U

1070 mm x 100 mm

Single

40270-U

TOURNIQUET CUFF COVERS

Tourniquet Covers

Small (box of 30)

40300

Tourniquet Covers

Large (box of 30)

40310

FOAM CUFF LINERS

Packet of 10

40320

EXSANGUINATORS

Rhys Davies

40400

Esmarch Bandages (5.5m)

L1112-55

SAFETY CONNECTORS

Male/Female complete

40600

Female only

40610

Male only

40620

Plastic Strip Ties

(pack of 50)

40630

O' Ring Washers

(pack of 50)

40640

Y CONNECTORS

For Compressed Air 4 Bar

40500

For Compressed Air 7 Bar

40510